

# IOL-Injector Correlation Chart

## IOL View



Seeing beyond

| IOL           |                    | Diopter Range & Injector / Cartridge                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|---------------|--------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4-Haptic MICS | CT SPHERIS® 209M   |    | <b>0.0 to +32.0 D</b><br>AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set<br>or VISCOJECT™-BIO 1.8 Injector Set<br>or A6 injector Set                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|               | CT ASPHINA® 409M   |    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|               | CT ASPHINA® 509M   |    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|               | AT LISA® 809M      |    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|               | AT TORBI® 709M     |  | <b>-10.0 to +28.0 D sphere; +1.0 to +4.0 D cylinder</b><br>AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set<br>or VISCOJECT™-BIO 1.8 Injector Set<br>or A6 injector Set<br><br><b>-10.0 to +28.0 D sphere; +4.5 to +12.0 D cylinder</b><br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set<br><br><b>+28.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder</b><br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set                                                           |
|               | AT TORBI® 719M     |  | <b>-9.5 to +28.0 D spherical equivalent (SE)*; +1.0 to +4.0 D cylinder</b><br>AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set<br>or VISCOJECT™-BIO 1.8 Injector Set<br>or A6 injector Set<br><br><b>-7.5 to +28.0 D spherical equivalent (SE)*; +4.5 to +12.0 D cylinder</b><br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set<br><br><b>+28.5 to +38.0 D spherical equivalent (SE)*; +1.0 to +12.0 D cylinder</b><br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set |
|               | AT LISA toric 909M |  | <b>-10.0 to +24.0 D sphere; +1.0 to +12.0 D cylinder</b><br>AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set<br>or VISCOJECT™-BIO 1.8 Injector Set<br>or A6 injector Set<br><br><b>+24.5 to +32.0 D sphere; +1.0 to +12.0 D cylinder</b><br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set                                                                                                                                                                                                                                  |




\*Further preselected SE/cylinder combinations are available above and below the stated SE range.

# IOL-Injector Correlation Chart

IOL View



Seeing beyond

| IOL           |                                 | Diopter Range & Injector / Cartridge                                                |                                                                                                                                                                                                                                                                       |
|---------------|---------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4-Haptic MICS | AT LISA tri toric® 939M         |    | <b>-10.0 to +28.0 D sphere;</b> +4.5 to +12.0 D cylinder and<br><b>+28.5 to +32.0 D sphere;</b> +1.0 to +12.0 D cylinder<br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set                                      |
|               | AT LISA tri toric® 949M         |    | <b>-5.0 to +32.0 D spherical equivalent (SE);</b> +4.5 to +12.0 D cylinder and<br><b>+32.5 to +35.0 D spherical equivalent (SE);</b> +1.0 to +12.0 D cylinder<br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set |
|               | AT LARA toric® 929M             |    | <b>-4.0 to +32.0 D spherical equivalent (SE)*;</b> +4.5 to +12.0 D cylinder<br>VISCOJECT™-BIO 2.2 Injector Set<br>or AT.Shooter A1-2000 injector / VISCOJECT™-BIO 2.2 Cartridge Set                                                                                   |
|               | <b>BLUEMIXS® Preloaded IOLs</b> |                                                                                     |                                                                                                                                                                                                                                                                       |
|               | CT ASPHINA® 409MP               |  | <b>0.0 to +32.0 D</b><br>BLUEMIXS® 180 Injector                                                                                                                                                                                                                       |
|               | CT ASPHINA® 419MP               |  |                                                                                                                                                                                                                                                                       |
|               | CT ASPHINA® 509MP               |  |                                                                                                                                                                                                                                                                       |
|               | AT LISA 809MP                   |  |                                                                                                                                                                                                                                                                       |
|               | AT LISA tri® 839MP              |  |                                                                                                                                                                                                                                                                       |
|               | AT LARA® 829MP                  |  | <b>-10.0 to +28.0 D sphere;</b> +1.0 to +4.0 D cylinder<br>BLUEMIXS® 180 Injector                                                                                                                                                                                     |
|               | AT TORBI® 709MP                 |  |                                                                                                                                                                                                                                                                       |
|               | AT TORBI® 719MP                 |  | <b>-8.0 to +28.0 D spherical equivalent (SE)*;</b> +1.0 to +4.0 D cylinder<br>BLUEMIXS® 180 Injector                                                                                                                                                                  |
|               | AT LISA tri toric® 939MP        |  | <b>-10.0 to +28.0 D sphere;</b> +1.0 to +4.0 D cylinder<br>BLUEMIXS® 180 Injector                                                                                                                                                                                     |
|               | AT LISA tri toric® 949MP        |  | <b>-5.0 to +32.0 D spherical equivalent;</b> +1.0 to +4.0 D cylinder<br>BLUEMIXS® 180 Injector                                                                                                                                                                        |

\*Further preselected SE/cylinder combinations are available above and below the stated SE range.

# IOL-Injector Correlation Chart

## IOL View



Seeing beyond

| IOL      |                      | Diopter Range & Injector / Cartridge                                              |                                                                                                                        |
|----------|----------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| 3-Haptic | AT LARA® toric 929MP |  | <b>-8.0 D to +32.0 D spherical equivalent (SE)*</b> ; +1.0 to +4.0 D cylinder<br>BLUEMIXS® 180 Injector                |
|          | CT SPHERIS® 204      |  | <b>-10.0 to +30.0 D</b><br>AT.Shooter A1-2000 / VISCOJECT™-BIO 2.2 Cartridge Set<br>or VISCOJECT™-BIO 2.2 Injector Set |
|          | CT ASPHINA® 404      |  | <b>+31.0 to +45.0 D (CT SPHERIS 204)</b><br><b>+31.0 to +42.0 D (CT ASPHINA 404)</b><br>ACCUJECT™ 3.0 Injector Set     |

\*Further preselected SE/cylinder combinations are available above and below the stated SE range.

# IOL-Injector Correlation Chart

## IOL View



Seeing beyond

| IOL            |                            | Diopter Range & Injector / Cartridge                                                |                                                                                                                                                                               |
|----------------|----------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1-Piece C-Loop | Preloaded hydrophobic IOLs |                                                                                     |                                                                                                                                                                               |
|                | CT LUCIA® 621P/PY          |    | <b>0.0 to +24.0 D</b><br>BLUESERT™ 2.2 Injector Set<br><br><b>+24.5 to +30.0 D</b><br>BLUESERT™ 2.4 Injector Set<br><br><b>+30.5 to +34.0 D</b><br>BLUESERT™ 2.6 Injector Set |
|                | CT LUCIA® 221P             |    | <b>0.0 to +24.0 D</b><br>BLUESERT™ 2.2 Injector Set<br><br><b>+24.5 to +30.0 D</b><br>BLUESERT™ 2.4 Injector Set<br><br><b>+30.5 to +34.0 D</b><br>BLUESERT™ 2.6 Injector Set |
|                | AT ELANA® 841P             |  | <b>0.0 to +24.0 D</b><br>BLUESERT™ 2.2 Injector Set<br><br><b>+24.5 to +30.0 D</b><br>BLUESERT™ 2.4 Injector Set<br><br><b>+30.5 to +34.0 D</b><br>BLUESERT™ 2.6 Injector Set |

# IOL-Injector Correlation Chart

## IOL View



Seeing beyond

| IOL     |                                              | Diopter Range & Injector / Cartridge                                              |                                                                                      |
|---------|----------------------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 3 Piece | CT LUCIA® 202                                |  | <b>+4.0 to +27.0 D</b><br>ACCUJECT™ 3.0 3P<br><br><b>+27.5 to +34.0 D</b><br>Forceps |
|         | CT LUCIA® 602<br>(Only available in the USA) |  | <b>+4.0 to +34.0 D</b><br>R28 Injector / Z28 Cartridge                               |

# IOL-Injector Correlation Chart

## Injector View



Seeing beyond


| Injector                                       | Injector                                           |                                                                                      | Injector Type            | Incision Size |
|------------------------------------------------|----------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------|---------------|
| AT.Shooter<br>A1-2000                          | VISCOJECT™-BIO 1.8<br>Cartridge Set <sup>1,2</sup> |    | Reusable                 | 1.8 mm        |
| VISCOJECT™-BIO 1.8 Injector Set <sup>2,3</sup> |                                                    |    | Disposable               |               |
| Single Use A6 Injector Set                     |                                                    |   | Disposable               |               |
| BLUEMIXS® 180                                  |                                                    |  | Disposable,<br>Preloaded |               |
| AT.Shooter<br>A1-2000                          | VISCOJECT™-BIO 2.2<br>Cartridge Set <sup>1,2</sup> |  | Reusable                 | ≥ 2.2 mm      |
| VISCOJECT™-BIO 2.2 Injector Set <sup>2</sup>   |                                                    |  | Disposable               |               |

# IOL-Injector Correlation Chart

## Injector View



Seeing beyond

| Injector                   | Injector                   |                                                                                                                                                                              | Injector Type | Incision Size |
|----------------------------|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|---------------|
| BLUESERT™ 2.2 <sup>5</sup> |                            |                                                                                            | Disposable    | 2.2 mm        |
| BLUESERT™ 2.4              |                            |                                                                                                                                                                              | Disposable    | 2.4 mm        |
| BLUESERT™ 2.6              |                            |                                                                                                                                                                              | Disposable    | 2.6 mm        |
| R28 Injector <sup>4</sup>  | Z28 Cartridge <sup>4</sup> |                                                                                            | Reusable      | ≥ 2.8 mm      |
| ACCUJECT™ 3.0 1P           |                            | <br> | Disposable    | ≥ 3.0 mm      |
| ACCUJECT™ 3.0 3P           |                            |                                                                                                                                                                              |               |               |

1 Soft silicone cushion is included in the cartridge set

2 VISCOJECT™-BIO or "classic" versions

3 VISCOJECT™ 1.8 ≥ 1.8 mm with docking technique

4 Only available in the USA

5 Not sold as stand alone